

Biotinylated rAB Expression



1.0 Introduction

This method is for the mid-scale expression and purification of site-specifically biotinylated Fabs from the RH2.2-Avi pTAC-based IPTG-inducible, Avi-tagged vector.

For the most up to date version of this protocol please visit: <http://recombinant-antibodies.org/protocols/rh22-avi>

2.0 Materials

Glassware/Plasticware

- 250 mL baffled flasks (Fisher Scientific, 10-140-6A)
- 50 mL Conical Tube (Bioexpress, LLC C-3394-4)
- 15 mL Biorad protein purification columns (EconoPac, 732-1010)
- 1.7 mL Eppendorf tubes (Axygen, MCT-150-C)

Reagents

Name	Company	Catalog Number
N-Z Amine	Fluka	C02626-1KG
yeast extract	Bioshop	YEX401.1
bio-tryptone	Bioshop	TRP402.1
glucose	Sigma	G8270-1KG
carbenicillin (carb)	Bioshop	CAR544.10
chloramphenicol (cmp)	Bioshop	CLR201.100
biotin	Sigma	B4501-5G
Tween 20	Bioshop	TWN510.500
potassium chloride	Bioshop	POC308.1
calcium chloride	Bioshop	CCL302.500
magnesium chloride	Bioshop	MAG510.500
monobasic potassium phosphate	Bioshop	PPM302.1
dibasic sodium phosphate	Anachemia	84486-440
ammonium chloride	Amresco	0621-1KG
ammonium sulphate	Bioshop	AMP302
sodium sulphate	Bioshop	SOS513.500
magnesium sulphate	BDH	BDH0246-500G

sodium chloride	Bioshop	Sod001.10
phosphoric acid	Acros Organics	201140010
lysozyme	Bioshop	LYS702.25
benzonase	Novagen	71205
Triton X-100	Bioshop	TRX506.500
protease inhibitor cocktail tablets	Roche	11 836 170 001
Protein A sepharose	GE Healthsciences	17-1279-03
IPTG	Bioshop	IPT001.50

BirA-transformed, chemically competent T1^r BL21 *E. coli*

Solutions

1X KCM

KCl (500 mM)

CaCl₂ (150 mM)

MgCl₂ (250 mM)

2YT media

1% yeast extract (10 g)

1.6% bio-tryptone (16 g)

86 mM NaCl (5 g)

Make up volume to 1L and adjust pH to 7.0

20X salt stock

500 mM Na₂HPO₄ (70.98 g)

500 mM KH₂PO₄ (68.05 g)

1M NH₄Cl (53.4 g)

100 mM Na₂SO₄ (14.20 g)

Mix with heating and sterile filter or autoclave until dissolved and make up to 1L. Store at RT

NZY media

1% N-Z-Amine (10g)

0.5% yeast extract (5 g)

Mix in 500 mL and stir until dissolved

Add 50 mL of 20X salt

Add 2 mL of 1M MgSO₄

Make up volume to 980 mL and autoclave

Lysis buffer:

50mM Tris

150mM NaCl

1% Triton X-100

1mg/ml lysozyme

2mM MgCl₂

10U benzonase

Elution buffer:

50mM NaH₂PO₄

140mM NaCl

100mM H₃PO₄

pH 2.5

3.0 Transformation

3.1 Thaw 25 µL of chemically competent BL21-BirA T1^r cells on ice

3.2 Add 50 ng of sequence-verified pTAC-based RH2.2_Avi-expression plasmid to a mixture of 4 µL of 5X KCM in 16 µL of MQ H₂O on ice

3.3 Chill mixture on ice for 10 minutes

3.4 Add 20 mL of chemically competent BL21-BirA T1^r cells to the mixture of KCM and DNA

3.5 Incubate 20 min on ice, transfer to the benchtop, incubate 10 min at RT, return to ice and incubate 2 min on ice

4.0 Expression

4.1 Following the last stage of transformation on ice, transfer the entirety of the cell/DNA mixture in to 25 mL of 2YT supplemented with 50 µg/mL carb and 5 µg/mL cmp in a 50 mL Falcon tube, leaving the lid loose but taped secure to allow gas transfer

4.2 Inoculate 1/40th volume (2.5 mL) of overnight culture into NZY media supplemented with 50 µg/mL carb + 5 µg/mL cmp + 25 µM biotin (100 µL) and 1/20th volume of 20X salt mixture

4.3 Grow 2-3 hr at 37°C, 200rpm until OD_{0.8-1.0} achieved

4.4 Induce culture with 1/1000th volume of 1M IPTG (i.e. 100 µL in to 100 mL of culture, final [IPTG] = 1mM)

4.5 Continue growing for 6-8 hrs at 30°C and shaking at 200 rpm

4.6 Transfer ½ the volume of culture to a 50 mL Falcon tube and pellet cells at 8000 rpms for 20 min.

4.7 Discard supernatant and repeat to obtain a single pellet.

4.6 Discard remaining supernatant. At this point the pellet can be frozen for purification at a later date or lysed for purification as per the following step.

5.0 Purification

5.1 Lyse pellet by the addition of 10 mL of lysis buffer and nutation for 1.5-4 hrs at 4°C

5.2 To ensure complete lysis, mixture can be sonicated for 1 min at 40% intensity using a 1 min program of 5 seconds on and 5 seconds off.

5.2 Centrifuge the pellet lysate at 9000rpm for 20min.

5.3 Add 250 µl of a protein A sepharose slurry added to a 15 ml GF column.

5.4 Equilibrate resin with 25 mL of 1X PBS, drain and cap.

5.5 Transfer the lysate supernatant to the equilibrated protein A column and allow binding to proceed for 30 min at 4°C with occasional mixing of the resin with the supernatant.

5.6 Drain the supernatant from the column retaining the flow-through and wash the resin bed with 40 ml of 1X PBS.

5.7 Add 60 µl of 1M Tris pH 8 added to Eppendorf collection tubes to neutralize the elution buffer.

4.8 Add 300 μ l Fab elution buffer (above) to each column and incubate 5 min before draining in to 1.7 mL Eppendorf containing neutralization buffer.

5.9 Elution step can be repeated for a total elution volume of 600 μ l neutralized with 1200 μ l Tris.

5.10 Additional polishing steps can be carried out as necessary.

5.11 Columns can be regenerated with 10mL of 100mM H_3PO_4 , washed with 20ml TBS, stored in 2ml 20% EtOH at 4°C.

Please send corrections, modifications and suggestions to shane.miersch@utoronto.ca

Contact:

Recombinant Antibody Network

admin@recombinant-antibodies.org

The University of Chicago

Knapp Center for Biomedical Discovery Rm. 3240G

900 E. 57th St., Chicago, IL 60637

Phone: +1 (773) 834-2776

University of California, San Francisco

Byers Hall Rm. 503

1700 4th St., San Francisco, CA 94158

Phone: +1 (530) 341-2371

University of Toronto

Best Institute Rm. 117

112 College Avenue, Toronto, Ontario M5G 1L6

Phone: +1 (416) 978-1594