

## 1.0 Introduction

This method is for the mid-scale expression and purification of Fabs from the RH2.2 pTAC-based IPTG inducible vector.

For the most up to date version of this protocol please visit: <http://recombinant-antibodies.org/protocols/rh22>

## 2.0 Materials

### Glassware/Plasticware

250 mL baffled flasks (Fisher Scientific, 10-140-6A)

50 mL Conical Tube (Bioexpress, LLC C-3394-4)

15 mL Biorad protein purification columns (EconoPac, 732-1010)

1.7 mL Eppendorf tubes (Axygen, MCT-150-C)

### Reagents

Name	Company	Catalog Number
N-Z Amine	Fluka	C02626-1KG
yeast extract	Bioshop	YEX401.1
bio-tryptone	Bioshop	TRP402.1
glucose	Sigma	G8270-1KG
carbenicillin	Bioshop	CAR544.10
Tween 20	Bioshop	TWN510.500
potassium chloride	Bioshop	POC308.1
calcium chloride	Bioshop	CCL302.500
magnesium chloride	Bioshop	MAG510.500
monobasic potassium phosphate	Bioshop	PPM302.1
dibasic sodium phosphate	Anachemia	84486-440
ammonium chloride	Amresco	0621-1KG
ammonium sulphate	Bioshop	AMP302
sodium sulphate	Bioshop	SOS513.500
magnesium sulphate	BDH	BDH0246-500G
sodium chloride	Bioshop	Sod001.10
phosphoric acid	Acros Organics	201140010
lysozyme	Bioshop	LYS702.25

benzonase	Novagen	71205
Triton X-100	Bioshop	TRX506.500
protease inhibitor cocktail tablets	Roche	11 836 170 001
Protein A sepharose	GE Healthsciences	17-1279-03
IPTG	Bioshop	IPT001.50
chemically competent T1 <sup>r</sup> BL21 <i>E. coli</i>		

## Solutions

### 1X KCM

KCl (500 mM)

CaCl<sub>2</sub> (150 mM)

MgCl<sub>2</sub> (250 mM)

### 2YT media

1% yeast extract (10 g)

1.6% bio-tryptone (16 g)

86 mM NaCl (5 g)

Make up volume to 1L and adjust pH to 7.0

### 20X salt stock

500 mM Na<sub>2</sub>HPO<sub>4</sub> (70.98 g)

500 mM KH<sub>2</sub>PO<sub>4</sub> (68.05 g)

1M NH<sub>4</sub>Cl (53.4 g)

100 mM Na<sub>2</sub>SO<sub>4</sub> (14.20 g)

Mix with heating and sterile filter or autoclave until dissolved and make up to 1L. Store at RT

### NZY media

1% N-Z-Amine (10g)

0.5% yeast extract (5 g)

Mix in 500 mL and stir until dissolved

Add 50 mL of 20X salt

Add 2 mL of 1M MgSO<sub>4</sub>

Make up volume to 980 mL and autoclave

### Lysis buffer:

50mM Tris

150mM NaCl

1% Triton X-100

1mg/ml lysozyme

2mM MgCl<sub>2</sub>

10U benzonase

#### Elution buffer:

50mM NaH<sub>2</sub>PO<sub>4</sub>

140mM NaCl

100mM H<sub>3</sub>PO<sub>4</sub>

pH 2.5

## 3.0 Transformation

3.1 Thaw 25 µL of chemically competent BL21-T1<sup>r</sup> cells on ice

3.2 Add 50 ng of sequence-verified pTAC-based RH2.2 expression plasmid to a mixture of 4 µL of 5X KCM in 16 µL of MilliQ H<sub>2</sub>O on ice

3.3 Chill mixture on ice for 10 minutes

3.4 Add 20 mL of chemically competent BL21-T1<sup>r</sup> cells to the mixture of KCM and DNA

3.5 Incubate 20 min on ice, transfer to the benchtop, incubate 10 min at RT, return to ice and incubate 2 min on ice

## 4.0 Expression

4.1 Following the last stage of transformation on ice, transfer the entirety of the cell/DNA mixture in to 25 mL of 2YT supplemented with 50 µg/mL carb in a 50 mL Falcon tube, leaving the lid loose but taped secure to allow gas transfer

4.2 Inoculate 1/40<sup>th</sup> volume (2.5 mL) of overnight culture into NZY media + 50 µg/mL carb (100 mL) supplemented with 1/20<sup>th</sup> volume of 20X salt mixture

4.3 Grow 2-3 hr at 37°C, 200rpm until OD 0.8-1.0 achieved

4.4 Induce culture with 1/1000<sup>th</sup> volume of 1M IPTG (i.e. 100 µL in to 100 mL of culture)

4.5 Continue growing for 6-8 hrs at 30°C and shaking at 200 rpm

4.6 Transfer ½ the volume of culture to a 50 mL Falcon tube and pellet cells at 8000 rpms for 20 min.

4.7 Discard supernatant and repeat to obtain a single pellet.

4.8 Discard remaining supernatant. At this point the pellet can be frozen for purification at a later date or lysed for purification as per the following step.

## 5.0 Purification

5.1 Lyse pellet by the addition of 10 mL of lysis buffer and nutation for 1.5-4 hrs at 4°C

5.2 To ensure complete lysis, mixture can be sonicated for 1 min at 40% intensity using a 1 min program of 5 seconds on and 5 seconds off.

5.3 Centrifuge the pellet lysate at 9000rpm for 20min.

5.4 Add 250 µl of a protein A sepharose slurry to a 15 ml GF column.

5.5 Equilibrate resin with 25 mL of 1X PBS, drain and cap.

5.6 Transfer the lysate supernatant to the equilibrated protein A column and allow binding to proceed for 30 min at 4°C with occasional mixing of the resin with the supernatant.

5.7 Drain the supernatant from the column retaining the flow-through and wash the resin bed with 40 ml of 1X PBS.

5.8 Add 60 µl of 1M Tris pH 8 added to Eppendorf collection tubes to neutralize the elution buffer.

5.9 Add 300 µl Fab elution buffer (above) to each capped column and incubate 5 min before draining in to an Eppendorf containing neutralization buffer.

5.10 Elution step can be repeated for a total elution volume of 600 µl neutralized with 1200 µl Tris.

5.11 Additional polishing steps can be carried out as necessary.

5.12 Columns can be regenerated with 10mL of 100mM  $\text{H}_3\text{PO}_4$ , washed with 20ml TBS, stored in 2ml 20% EtOH at 4°C.

Please send corrections, modifications and suggestions to [shane.miersch@utoronto.ca](mailto:shane.miersch@utoronto.ca)

**Contact:**

**Recombinant Antibody Network**

admin@recombinant-antibodies.org

**The University of Chicago**

Knapp Center for Biomedical Discovery Rm. 3240G

900 E. 57th St., Chicago, IL 60637

Phone: +1 (773) 834-2776

**University of California, San Francisco**

Byers Hall Rm. 503

1700 4th St., San Francisco, CA 94158

Phone: +1 (530) 341-2371

**University of Toronto**

Best Institute Rm. 117

112 College Avenue, Toronto, Ontario M5G 1L6

Phone: +1 (416) 978-1594